

PROTECTION SET

SET FOR OPERATOR, ASSISTANT AND PATIENT

3 different level of protection, 3 different material quality

The individually packed sets are composed in configuration ideal for achieving a good level of isolation. Each set includes the necessary disposables to protect operator, assistant and patient during the daily procedures.



10.T0799.00 for operator*
10.T0798.00 for assistant*
Premium High Protection Set



10.T0564.00 for operator*
10.T0594.00 for assistant*
High Protection Set



10.T0600.00 for operator*
10.T0601.00 for assistant*
Medium Protection Set



10.S9378.00*
Patient Protection Set

DESCRIPTION

Set operator and assistant

Each Set has been designed to dress the operator and the assistant. It contains a mask a gown and a cap to ensure a perfect isolation during your entire day, in order to get a high level of protection to better compensate the environment deficiencies. Just wear it for your first patient and then, in combination with the gown cover present into the patient Set, change just it for every new patient.

Set patient

Each set has been designed to dress and protect the patient. It also contains 2 gown covers for operator and assistant to wear every new patient treatment.

Daily use: 1 each patient

*Available also sterile.

FEATURES AND BENEFITS

- **Compliant** to the stricter national and local regulations.
- More attention on operators' and staff's **safety**.
- Increased **patients'** attention.
- **Timesaving**, better **organization** and **efficiency** for the chairside preparation.

PREMIUM HIGH PROTECTION SET

For operator and assistant

Premium High Protection Set - Operator ref. 10.T0799.00 10 pz/cf			Premium High Protection Set - Assistant ref. 10.T0798.00 10 pz/cf		
Description	Ref. Code	Qty	Description	Ref. Code	Qty
Isolator™ N95 Particulate Respirator	GPRN95	1	Isolator™ N95 Particulate Respirator	GPRN95	1
SMS/PP Water resistant helmet-like cap with laces	30.M0046.00	1	SMS/PP Water resistant helmet-like cap with laces	30.M0046.00	1
Surgically folded model gown - XL wrapped in medical paper with 2 hand towels	30.D1276.00	1	Surgically folded model gown - S wrapped in medical paper with 2 hand towels	30.D1278.00	1
Face shield	10.M0200.00	1	Face shield	10.M0200.00	1

Available also sterile: ref. 12.T0799.00

Available also sterile: ref. 12.T0798.00



SMS/PP Water resistant helmet-like cap with laces

Water / Fluid resistant - Breathable part.



Surgically folded model gown, wrapped in medical paper with 2 hand towels

- Spunlace medical nonwoven fabric.
- High liquid-repellence.
- Extraordinary sense of comfort.
- Velcro closure, cotton cuffs.
- Double internal and external side laces.
- Surgically folded.
- 2 hand towels.
- Wrapped in medical paper.



Isolator™ N95 Particulate Respirator

- NIOSH approved N95 Respirator.
- Meets guidelines for TB exposure control.
- Resistance to penetration by synthetic blood: 160 mmHg.
- Inner/outer fluid resistant layers.
- Chin pocket construction.
- EN 149: 2001 FFP2 NR / NIOSH N95

Although the N95 closely resembles a surgical mask, it is a respirator and must be used in accordance with all local regulations regarding respiratory protection.



Face shield

- Wide polyethylene fog preventing shield.
- Polyurethane insert.
- Suitable by people wearing glasses.

**THE ONLY DIFFERENCE
BETWEEN THE OPERATOR
SET AND THE ASSISTANT SET
IS THE SIZE OF THE GOWN**

HIGH PROTECTION SET

For operator and assistant

High Protection Set - Operator - ref. 10.T0564.00 - 10 pz/cf

Description	Ref. Code	Qty
Isolator™ N95 Particulate Respirator	GPRN95	1
SMS Water resistant helmet-like cap with laces	30.M0046.00	1
Surgically folded gown - wrapped in medical paper (length 118 cm) XL with 2 hand towels	30.D1273.00	1

Available also sterile: ref. 12.T0564.00



Isolator™ N95 Particulate Respirator

- NIOSH approved N95 Respirator.
- Meets guidelines for TB exposure control.
- Resistance to penetration by synthetic blood: 160 mmHg.
- Inner/outer fluid resistant layers.
- Chin pocket construction.
- EN 149: 2001 FFP2 NR / NIOSH N95

Although the N95 closely resembles a surgical mask, it is a respirator and must be used in accordance with all local regulations regarding respiratory protection.

Surgically folded model gown, wrapped in medical paper with 2 hand towels

- Fluid-repellent and breathable SMMS Microfibre.
- Excellent comfort/performance ratio.
- Velcro closure - cotton cuffs.
- Side-fastened waist-belt.
- Surgical folding.
- Internal 2 towels.
- Wrapped in medical paper.



High Protection Set - Assistant - ref. 10.T0594.00 - 10 pz/cf

Description	Ref. Code	Qty
Isolator™ N95 Particulate Respirator	GPRN95	1
SMS/PP Water resistant helmet-like cap with laces	30.M0046.00	1
Surgically folded gown - wrapped in medical paper (length 112 cm) M with 2 hand towels	30.D1271.00	1

Available also sterile: ref. 12.T0594.00



SMS/PP Water resistant helmet-like cap with laces

Water / Fluid resistant - Breathable part.

MEDIUM PROTECTION SET

For operator and assistant

Medium Protection Set - Operator ref. 10.T0600.00 - 10 pz/cf

Description	Ref. Code	Qty
Ultra™ Sensitive Earloop Mask with secure Fit™ Mask Technology	GCFCXSFSF	1
Surgical Cap with adjustable lace	30.M0045.00	1
Surgically folded SMMS model gown (35 gr) size L	30.D2302.00	1

Available also sterile: ref. 12.T0600.00



Ultra™ Sensitive Earloop Mask with secure Fit™ Mask Technology

- Type IIR Enhanced Fluid Resistance.
- ASTM 3.
- Fluid resistant outer layer.
- Fluid resistant, extra soft white hypoallergenic inner cellulose layer – will not lint, tear or shred.
- Free of chemicals, inks and dyes.
- Extra-long aluminum nose piece.
- FogFree™ Strip on mask inner layer.

Surgical Cap with adjustable lace

- SMS fluid repellent.
- Enables ventilation.



Medium Protection Set - Assistant ref. 10.T0601.00 - 10 pz/cf

Description	Ref. Code	Qty
Ultra™ Sensitive Earloop Mask with secure Fit™ Mask Technology	GCFCXSFSF	1
Surgical Cap with adjustable lace	30.M0045.00	1
Surgically folded SMMS model gown (35 gr) size L	30.D2302.00	1

Available also sterile: ref. 12.T0601.00



Surgically folded SMMS model gown (35 gr)

- Water-repellent and breathable, SMMS.
- Cotton cuffs.
- Waist-belt.
- Surgical folding make.
- Ultrasonic welded seams to guarantee a more efficient barrier against bacteria and fluids.

PATIENT PROTECTION SET

For patient

Patient protection Set - ref. 10.S9378.00 - 10 pz/cf

Description	Ref. Code	Qty
Gown Cover - LDPE	10.M1002.00	2
Patient Safety Mask with shield	GCPAT	1
Patient cap	30.M0040.00	1
Sani-Tab™ Chain-Free Patient Towel	WEXBLT	1
Patient Shoe Cover	30.D4003.00	2
Round 24 cm Ø waterproof headrest with elastic band - white	20.O0011.00	1
Comfort Plus™ Premium Saliva Ejectors White/Green Mint scented	ZWGCPM	1
Adhesive films 20x20 cm LDPE - BLUE	80.F0008.00	2

Available also sterile: ref. 12.S9378.00

The two gown covers are intended to be used by operator and assistant during **every new patient treatment**.

DAILY USE:
1 per patient




Gown Cover - LDPE

- Medical grade
- Waterproof
- High insulation standard



Patient cap

- Fluid repellent.
- Enables ventilation.
- Regarding respiratory protection.



Sani-Tab™ Chain-Free Patient Towel

- Easy to use – simply position towel anywhere on the patient.
- Bib stays where placed to maximize coverage, unlike bib clips and other.
- disposable towel holders.
- Can be placed high on patient to protect clothing.




Patient Shoe Cover

Waterproof shoe s with elastic opening.



Round 24 cm Ø waterproof headrest with elastic band - white

Waterproof and absorbent NWF protections.



Comfort Plus™ Premium Saliva Ejectors White/Green Mint scented

- The kinder, gentler saliva ejector..takes in fluid, not tissue!
- Patented tip design enhances patient comfort – will not grab mucosal tissue.
- Facilitates uninterrupted fluid and light debris removal.
- Optimizes fluid removal as patented tip design prevents tissue from blocking suction.
- Holds shape when formed as result of proprietary manufacturing process.
- Respiratory protection.



Adhesive films 20x20 cm LDPE - BLUE

- LDPE protections.
- Avoid cross-contamination.



Patient Safety Mask with shield

- Fluid resistant. Covers patient's eyes and nose.
- Protects fragile eye and mucosal tissue from flying particles, spray and splatter.
- Makes patients feel safe, but not "in the dark," with attached clear-vision shield.
- Form-fitting profile, flexible materials provide maximum access to mouth.
- White spunbond inner and outer layers.

Manufacturer: Omnia srl | Via F. Delnevo 190 | 43036 Fidenza (PR) | Italy
Customer Care Department: Omnia srl
Tel +39 (0524 527453) | **Fax +39** (0524 525230)
E-Mail: info@omniaspa.eu | **Website:** omniaspa.eu

Hu-Friedy, Omnia, and Crosstex are now proud members of

HuFriedyGroup
The Best In Practice