

HuFriedyGroup

Key Opinion Leader

Prof. Francesco Cairo

Advanced Mucogengival and Regenerative Surgery Kit
& Periodontal Surgery Basic Kit



Focus on: **Periodontal and Peri-Implant Plastic surgery**



Country: **Italy**



Favourite Instrument:
Orban knife KO1/2X



Website:
www.francescocairo.com



Inspired by:
Prof. Giovan Paolo Pini Prato



VISIT THE
KOL COMMUNITY




HuFriedyGroup instruments are of extraordinary quality, are incredibly handy and precise, and are indispensable for high-precision surgery.

Prof. Francesco Cairo



ADVANCED MUCOGENGIVAL AND REGENERATIVE SURGERY KIT

Periodontal and peri-implant plastic surgery techniques are very sensitive and difficult. The present kit was developed in collaboration with HuFriedyGroup for colleagues that intended to improve skills in these procedures, reducing surgical trauma and improving flaps and grafts stability. These elements are critical to obtain excellent and stable soft tissue at teeth and implants.



Visit us online at **HuFriedyGroup.eu**

ADVANCED MUCOGENGIVAL AND REGENERATIVE SURGERY KIT



DID YOU KNOW?

All hinged instruments – forceps, rongeurs, scissors, pliers, hemostats, orthodontic pliers, etc. – should be kept lubricated. Regular use of proper lubricants, like Instrument Lubricant Spray (**ILS**), will prevent rust, corrosion and stiff joints to ensure smooth operation. **All hinged instruments should be sterilised in the open position.**



SPNHDPVN
Swiss Perio Microsurgical
Needle Holder



SPSPV
Swiss Perio Microsurgical
Scissors



TPG2
Tissue Pliers Gerald
anatomical angled 18 cm



SPTPDAPV
Swiss Perio Microsurgical Plier,
anatomical

BLACK LINE

The Black Line is designed to optimise clinical outcomes by delivering efficiency throughout the entire perio and surgical procedure. The Black Line features a performance engineered coating for optimal edge retention, reduced light reflection afforded by the black matte finish and enhanced contrast and visual acuity at the surgical site and underlying tissue.



SPP24GSICX
Swiss Periosteal #SIC
Black Line



SPP24GSPX
Swiss Periosteal #24G
Black Line




SPPSCHLEEX
Swiss Perio Periosteal
#SCHL **Black Line**



CSEUWES6
Euwe Bone Scoop/Plugger
hdl #6

ADVANCED MUCOGINGIVAL AND REGENERATIVE SURGERY KIT



TKN1X
Tunneling Knife #1
Black Line

TKN2X
Tunneling Knife #2
Black Line

KO1/2X
Knife Orban #1/2 Orban,
Black Line

KPAX
Knife Allen End-Cutting
Intrasulcular, **Black Line**



MIR5DSHDE
High Definition Mirror #5
double-sided, European,
1pc./pack

MHE6
Mirror Handle hdl #6
European Thread



C36/376
Back Action Bone
Chisel #36/37 hdl #6



PCPUNC156
Probe UNC #UNC15 hdl #6
Qulix 1-2-3-...15



10-130-05E
Scalpel Handle #5E,
European handle



SHDPV
Scalpel Handle Steel
Microsurgical

ADVANCED MUCOGENGIVAL AND REGENERATIVE SURGERY KIT

15



SAS1/291E2
Curette Gracey Mini Five™ #1/2 hdl #9 oral/labial, grey

16



SAS7/89E2
Curette Gracey Mini Five™ #7/8 hdl #9 oral/buccal, purple

17



SAS11/1293E2
Curette Gracey Mini Five™ #11/12 hdl #9 mesial, orange

18



SAS13/1498E2
Curette Gracey Mini Five™ #13/14 hdl #9 distal, blue

CHOOSE THE COLOR.

Just add the color coding number to the last part of the partcode (replace " _ " with the number)

Red	●	0
Grey	●	1
White	●	2
Orange	●	3
Purple	●	4
Yellow	●	5
Lavender	●	6
Black	●	7
Blue	●	8
Green	●	9
Ocean Blue	●	11


WHY YOU SHOULD USE AN INSTRUMENT MANAGEMENT SYSTEM (IMS):

By standardising your instruments with IMS System, additional time and revenue is generated while extending the useful life of your Hu-Friedy instruments. Dental professionals state that by implementing IMS™ Cassettes, they save 5 to 10 minutes per procedure allowing the staff to spend more time on value-added activities*. In addition, safer transportation of soiled instruments reduces risk of injury and exposure to blood borne pathogens.



*When compared to single instrument decontamination. Based on market survey results.



IMEDIN18_S
IMS™ Infinity Series™ **DIN Surgical cassette for 18 instruments with an accessory area + 4 clips (2 x 11000) for hinged instruments.**
Dimensions: 290 x 184 x 34 mm



PERIODONTAL SURGERY BASIC KIT



Classical periodontal and implant surgery requires several specific steps in the execution of the single techniques, from the flap design, to the roots treatment and suture. The present kit was developed in collaboration with HuFriedyGroup for colleagues that routinely applied periodontal surgery and comprises instruments that easily support the periodontist.



PCPUNC156
Probe UNC #UNC15 hdl
#6 Qulix 1-2-3-...15



CRM2
Cheek Retractor Cawood



S14
Scissors LaGrange
#14 2-x angled 11,5 cm

PERIODONTAL SURGERY BASE KIT



PQ2N7
Birfurcation Probe Nabers #2N hdl #7 Qulix 3-6-9-12



PPR3S6
Periosteal Prichard #3 short hdl #6



PBQ6
Periosteal Benque hdl #6 Satin Steel™



NH5024CORT
Needle Holder Cortellini Micro Castro #5024 round handle, straight



DP18PS
Dressing Pliers #18 Perma Sharp™

PERMA SHARP™
INSTRUMENTS

Featuring tungsten carbide inserts bonded to surgical-grade stainless steel on both cutting edges for smooth, controlled cutting and maximum durability.



10-130-05E
Scalpel Handle #5E, European handle



KGF76
Knife Goldman-Fox #7 hdl #6 retromolar oval



KGF86
Knife Goldman-Fox #8 hdl #6 bukal-lingual




KO1/26
Knife Orban #1/2 hdl #6



TP42
42 Adson 1x2 Tissue Pliers, 12 cm

PERIODONTAL SURGERY BASE KIT




CO1
Bone Chisel Ochsenbein
#1, hdl #524

CO2
Bone Chisel Ochsenbein
#2, hdl #524

C36/376
Back Action Bone Chisel
#36/37 hdl #6


FS9/10S6
File Schluger #9S/10S
hdl #6



MIR5DSHDE
High Definition Mirror #5
double-sided, European,
1pc./pack

MHE6
Mirror Handle hdl #6
European Thread

S17
Scissors Iris #17
straight 11,5 cm




SG1/29LCE2
Curette Gracey
Labriola-Cortellini
#1/2 hdl #9

SG7/89LCE2
Curette Gracey
Labriola-Cortellini
#7/8 hdl #9

SG11/129LCE2
Curette Gracey
Labriola-Cortellini
#11/12 hdl #9

SG13/149LCE2
Curette Gracey
Labriola-Cortellini
#13/14 DE hdl #9

S333/3477
Scaler Kaplan
#333/334 hdl #7
Mini, green



IME14DIN5_
IMS™ Infinity Series™ 1/4DIN cassette for 5 instruments.
Dimensions: 72 x 184 x 34 mm



IMEDIN18_S
IMS™ Infinity Series™ DIN Surgical Cassette for 18
instruments with an accessory area + 4 clips (2 x 1000)
for hinged instruments.
Dimensions: 290 x 184 x 34 mm

PERIOCHALLENGE ACADEMY

Professor of Periodontology and Director Master Program in Periodontology and Implant Dentistry (accredited by EFP European Federation of Periodontology) at Florence University, Florence, Italy.

He is also Director of post Graduate Residency Program in Oral Surgery and Head of Research Unit in Periodontology and Periodontal Medicine in the same University. He is currently the President Elect of SIdP, Italian Society of Periodontology and Implantology. He authored several original articles on international peer-reviewed journals. Winner of important clinical research prizes, including the H. Goldman from SIdP, the Jaccard prize from EFP and R. Earl Robinson award from American Academy of Periodontology (AAP). He is member of the Journal of Clinical Periodontology editorial board and referee for the most important periodontal and implant journals.

www.francescocairo.com



COURSES:

- Annual course in Periodontology
- Advanced course in Periodontal and Peri-Implant Plastic Surgery
- International Masterclass on Periodontal and Peri-Implant Plastic Surgery

For more information or to place an order, call at 00800 48 37 43 39 or visit HuFriedyGroup.eu

Free Call: 00800 48 37 43 39 | **Free Fax:** 00800 48 37 43 40

E-Mail: info@hu-friedy.eu | **Website:** HuFriedyGroup.eu

HuFriedyGroupEurope

Hu-Friedy Mfg. Co., LLC • European Headquarters • Lyoner Str. 9 • 60528 Frankfurt am Main, Germany • HuFriedyGroup.eu
All company and product names are trademarks of Hu-Friedy Mfg. Co. LLC., its affiliates or related companies, unless otherwise noted.
©2023 Hu-Friedy Mfg. Co., LLC. All rights reserved. HFL-726/0223

HuFriedyGroup
The Best In Practice