

HuFriedyGroup

Key Opinion Leader

Dr. Daniele Cardaropoli

Periodontal Kit



Focus on:

Periodontal & Implant surgery



Country: **Italy**



Favourite Instrument: **70.H1519.00 CARDAROPOLI Bone-Compactator 15 mm/spoon**



Website: **www.proed.it**



Inspired by: **Prof. Myron Nevins**



VISIT THE
KOL COMMUNITY



In my surgical practice I want to work using only the best quality products

Dr. Daniele Cardaropoli




PERIODONTAL KIT

The complete kit provides a set of sharp surgical instruments which are essential during the periodontal conservative and resective-regenerative surgery.

This set can be used both for remodelling bone tissue and decontaminating periodontal pockets and for polishing root surfaces.

The millimeter probe, up to 15 mm, is essential for proper diagnosis and planning during the whole treatment.



Visit us online at **HuFriedyGroup.eu**

PERIODONTAL KIT



SM239E2
#23 Posterior Sickie Scaler,
#9



C36/376
36/37 Back-Action DE Chisel,
#6



SGF46
Goldman-Fox #4 Anterior
Scaler , #6



PCPUNC156
Color-Coded Probe #CP-15
UNC Metal SE #6 Handle



10-130-05E
Scalpel Handle
#5E European Style
Round, Matte Finish



CV6193
Carver Dis #6193, hdl #41



PPAELA
Periosteal, Allen Elevator,
Anterior , #6



KKN116
11 Kramer-Nevins Knife DE, #6



CM2/4S6
2/4 Molt DE Curette, #6, Short

PERIODONTAL KIT



10 **PH26M**
2 Hourigan Periosteal , Modified




11 **SPV**
Castroviejo
Microsurgical
Scissors 18 cm,
Standard



12 **NH5024R**
Perma Sharp™
Castroviejo Needle
Holder, Round 18 cm

PERMA SHARP™
INSTRUMENTS

Featuring tungsten carbide
inserts bonded to surgical-
grade stainless steel on
both cutting edges for
smooth, controlled cutting
and maximum durability.




13 **TPG1**
Gerald Tissue Pliers, Straight/Serrated, 18 cm



14 **70.H1519.00**
CARDAROPOLI
Bone-Compactor 15mm/spoon

INSTRUMENT DESIGNED
IN COLLABORATION WITH
DANIELE CARDAROPOLI

The instrument Carda # 4, finds its
application in the techniques of ridge
preservation for the filling of the socket
after extraction and in immediate post-
extractive implants for filling the gap
between bone and implant.



15 **S14SC**
14 LaGrange Scissors, Curved,
Super-Cut , 11.5 cm (4.5"), 1 serrated
blade, Standard

SUPER-CUT

Cutting blade honed/
serrated to „razor“
sharpness with opposing
blade finished with fine
serrations to avoid
tissue slippage.

CHOOSE THE COLOR

Just add the color coding
number to the last part of
the partcode (replace " _ "
with the number)

Red	●	0
Grey	●	1
White	●	2
Orange	●	3
Purple	●	4
Yellow	●	5
Lavender	●	6
Black	●	7
Blue	●	8
Green	●	9
Ocean Blue	●	11

WHY YOU SHOULD USE AN INSTRUMENT
MANAGEMENT SYSTEM (IMS):


By standardising your instruments with IMS System,
additional time and revenue is generated while extending
the useful life of your Hu-Friedy instruments. Dental
professionals state that by implementing IMS™ Cassettes,
they save 5 to 10 minutes per procedure allowing the
staff to spend more time on value-added activities*.
In addition, safer transportation of soiled instruments
reduces risk of injury and exposure to blood borne
pathogens.

*When compared to single instrument decontamination.
Based on market survey results.

DISCOVER
MORE AT IMS
WEBPAGE



IMEDIN18_S
IMS™ Infinity Series™ **DIN Surgical Cassette**
for 18 instruments, with an accessory
area + 4 clips (2 x I1000) for hinged instruments.
Dimensions: 290 mm x 184 mm x 34 mm



SUTURES

PTFE SUTURES
MONOFILAMENT
NON ABSORBABLE - WHITE

PTFE Surgical Suture (Polytetrafluoroethylene) are ideal for all implant, periodontal and bone grafting procedures where the use of a monofilament suture with low bacterial adhesion is desirable. Unlike other synthetic monofilament sutures, the material is exceptionally well tolerated in the oral cavity.



32.Z3161.00
6/0 Extra Reverse Cutting, 10 mm,
3/8 circle, 45 cm, 12 Pcs./box

“

The concept for success in periodontal and implant surgery is related to the ability of perfectly performing every single step of the procedure. Starting from this point, the flap design, the line of incision, the split or full-thickness elevation are not the only critical stages of the surgery, but also the management of the soft tissues during the final phases of the operation has to be performed at the best. For this reason, PTFE can be considered as the gold standard for the choice of the suturing material. The perfect tightness of the knot, the sliding, the biocompatibility, the little plaque retention represent basic characteristics during the choice of the suture, and the PTFE has all of them. Their ability to keep the flaps well closed and the incapability to keep the plaque are a plus in the post surgery period, so to limit any possible complication that may occur as the opening of the flaps with separation of the soft tissues or infection on the suture line. For these reasons, our personal preference is for the Omnia PTFE in all the operations of periodontal plastic and regene-rative surgery, ridge preservation and bone regenerative surgery.

Dr. Daniele Cardaropoli

”



32.Z3061.00
6/0 Extra Reverse Cutting, 12 mm,
3/8 circle, 45 cm, 12 Pcs./box

PGA SUTURES
MULTIFILAMENT - ABSORBABLE
SYNTHETICS - VIOLET

PGA sutures have synthetic origin and are made with a multi-thread consisting of braided Polyglycolic Acid, coated with a layer of Polycaprolactone copolymer. The woven structure offers a great ease of knotting and excellent hold of the knot also has a low bacterial adhesion. Suture resorption will be completed between the fiftieth and ninetieth day. The hold of the suture will be 65% of the total after 2 weeks and 35% of the total at the end of the third week. Thanks to the non-need to remove the stitches, it is a very versatile suture that lends itself to multiple uses in the field of oral surgery.

DANIELE CARDAROPOLI'S
OMNIA SET FOR IMPLANTOLOGY

“

We have been using Omnia procedure kits for several years both for periodontal and implant surgery. The possibility to have a direct confrontation with the company and the helpfulness of the people who work in it, allowed us to prepare custom sets, specifically designed for us, according to our demands and needs. In this way the preparation of the operative field has become a fast procedure, accurate and effective, perfectly executed by the staff of our practice. Sterility control is definitely the first step to be able to perform surgery successfully.

Dr. Daniele Cardaropoli

”



12.S7534.00
Dr. Daniele Cardaropoli's Omnia Set for Implantology

DESCRIPTION

1 x	Special Plus model gown (length 109 cm) size M
1 x	Special Plus model gown (length 118 cm) size L
1 x	Round cap with elastic
2 x	Envelope cap with Ominia logo
2 x	3-ply mask with round anallergic elastics light blue
2 x	Sheaths 120 x 7 cm with inserter Omnisleeve
1 x	Cannula, model Omniasurg
2 x	Light blue absorbent/waterproof drape 50 x 50 cm
1 x	Light blue absorbent/waterproof drape 50 x 75 cm
1 x	Absorbent/Waterproof drape 75x90 cm with adhesive U-shaped

PROED PROFESSIONAL EDUCATION

MASTERCLASS IN PERIODONTAL PLASTIC AND REGENERATIVE SURGERY

3-day Full Immersion program with slide presentation, video presentation, live-surgery and hands-on.


Where: Torino, Italy

Why: to improve the skills for the treatment of periodontal infra bony defects and gingival recessions

Info: www.proed.it or corsi@proed.it



VISIT THE
WEBSITE



Periodontology is the core of dentistry. Maintaining and restoring the health of tooth supporting tissues is the best way to ensure both their esthetic and functional integrity. This book systematically collects the most up-to-date knowledge on periodontal surgery.

DISCOVER
MORE



For more information or to place an order, call at 00800 48 37 43 39 or visit HuFriedyGroup.eu

Free Call: 00800 48 37 43 39 | **Free Fax:** 00800 48 37 43 40

E-Mail: info@hu-friedy.eu | **Website:** HuFriedyGroup.eu

f **in** **in** HuFriedyGroupEurope

Hu-Friedy Mfg. Co., LLC • European Headquarters • Lyoner Str. 9 • 60528 Frankfurt am Main, Germany • HuFriedyGroup.eu
All company and product names are trademarks of Hu-Friedy Mfg. Co. LLC., its affiliates or related companies, unless otherwise noted.
©2023 Hu-Friedy Mfg. Co., LLC. All rights reserved. HFL-782/0323

HuFriedyGroup
The Best In Practice